

REVIEW

Open Access

Supportive development of functional tissues for biomedical research using the MINUSHEET® perfusion system

Will W Minuth* and Lucia Denk

Abstract

Functional tissues generated under in vitro conditions are urgently needed in biomedical research. However, the engineering of tissues is rather difficult, since their development is influenced by numerous parameters. In consequence, a versatile culture system was developed to respond the unmet needs. Optimal adhesion for cells in this system is reached by the selection of individual biomaterials. To protect cells during handling and culture, the biomaterial is mounted onto a MINUSHEET® tissue carrier. While adherence of cells takes place in the static environment of a 24 well culture plate, generation of tissues is accomplished in one of several available perfusion culture containers. In the basic version a continuous flow of always fresh culture medium is provided to the developing tissue. In a gradient perfusion culture container epithelia are exposed to different fluids at the luminal and basal sides. Another special container with a transparent lid and base enables microscopic visualization of ongoing tissue development. A further container exhibits a flexible silicone lid to apply force onto the developing tissue thereby mimicking mechanical load that is required for developing connective and muscular tissue. Finally, stem/progenitor cells are kept at the interface of an artificial polyester interstitium within a perfusion culture container offering for example an optimal environment for the spatial development of renal tubules. The system presented here was evaluated by various research groups. As a result a variety of publications including most interesting applications were published. In the present paper these data were reviewed and analyzed. All of the results point out that the cell biological profile of engineered tissues can be strongly improved, when the introduced perfusion culture technique is applied in combination with specific biomaterials supporting primary adhesion of cells.

Keywords: Cell culture, Tissue culture, Tissue engineering, Biomaterial, MINUSHEET® tissue carrier, Perfusion culture, Gradient perfusion culture, Artificial polyester interstitium

Review

Background

In tissue engineering, biomaterial research and regenerative medicine culture experiments with epithelial, connective, muscular and nervous tissues are of extreme importance. In the focus of interest are cellular interactions with innovative biomaterials, controlled tissue regeneration and toxic influences of newly developed pharmaceuticals preventing inflammation but promoting renewal of tissues.

While the expansion of isolated cells in the static environment of a culture dish is easy to perform, developing tissues frequently illustrate severe morphological, physiological and biochemical alterations caused by dedifferentiation especially in combination with innovative biomaterials [1-8]. Atypical features are not evoked by one parameter, but are influenced on the one hand by suboptimal material surfaces resulting in minor cell adhesion and communication. On the other hand culture conditions such as suboptimal nutrition, low respiratory gas or unstirred layers of fluid are responsible for the process of atypical development [9-14]. Consequently, several parameters have to complement one another for an optimal generation of tissues. Since a classical dish

* Correspondence: will.minuth@vkl.uni-regensburg.de
Department of Molecular and Cellular Anatomy, University of Regensburg, Regensburg, Germany

does not comply to necessary needs of cells, advanced culture techniques allowing creation of an adaptable environment have to be applied for intact development of tissues [15].

Due to the multiple parameters that influence tissue development a modular system is presented here that adjusts culture conditions to individual needs of developing tissues [16]. Since suitable equipment for tissue engineering was at that time not commercially available, the described tools were developed for the own laboratory work. The innovative concept is based on a MINUSHEET® tissue carrier, which enables the researcher to apply individually selected biomaterials promoting in turn optimal cell adhesion. After transfer of a tissue carrier in a compatible perfusion culture container, the fluid environment can be individually macheted to the contained tissue. Thus, the modular system provides a highly adaptable basis for adhesion of cells in a MINUSHEET® tissue carrier with subsequent culture in the static environment of a 24 well culture plate. The compatible tissue carriers in combination with different kinds of perfusion culture containers made it possible to generate specialized tissues under continuous provision with always fresh medium.

Modulating tissue development by an adapted environment

Under natural conditions a basic prerequisite for optimal tissue development is a positive interaction of adherent cells with the extracellular matrix. Under in vitro conditions a selected biomaterial has to replace the naturally occurring extracellular matrix.

Balance between differentiation and dedifferentiation

Tissue engineering experiments demonstrated that the composition of applied materials such as synthetic polymers, decellularized matrix, biodegradable scaffolds, ceramics or metals has a strong influence on the adhesion of cells and consequently on the generation of tissues [17-26]. Furthermore, the roughness of an applied biomaterial can determine the development of maturing tissue in a good or bad sense. Also the spatial contour of biomaterials available in form of filters, nets, foils, fleeces or foams including small or big pores plays an essential role in regulation of tissue formation. Regarding the huge amount of newly developed biomaterials the challenge for the future is to recognize as early as possible the respective positive influences or toxic features on adherent cells by an advanced culture system.

Mounting biomaterials in a tissue carrier

Adherent cells can be cultured on individually selected biomaterials as a substitute for the extracellular matrix. To stay compatible with a 24 well culture plate, the

selected biomaterials are excised by a punching tool to a diameter of 13 mm or are ordered in this format (Figure 1a). For easy handling and to prevent damage during culture, the selected material is mounted in the base part of a MINUSHEET® tissue carrier (Figure 1b; black ring). The biomaterial is held in position by pressing it down with a tension ring (Figure 1c; white ring). The mounted tissue carrier is then enveloped and sterilized depending on the chemical composition of contained biomaterial either in formalin, ethylene oxide gas atmosphere or in an autoclave. Subsequently the tissue carrier can be frozen, stored at room temperature or used immediately for cell seeding.

Seeding of cells on a tissue carrier

For cell seeding the tissue carrier including the biomaterial of choice is placed by forceps in a 24 well culture plate (Figure 1d). To increase cell concentration, the selected biomaterial is solely wetted with culture medium. A standard culture on a MINUSHEET® tissue carrier can be performed by seeding cells onto the upper side of the selected biomaterial. When the carrier is turned, cells can be also seeded on the basal side so that co-culture with two different cell types becomes possible [27,28].

Cells are transferred to the selected biomaterial by a pipette within a small droplet of medium. Since the tissue carrier rests on small protrusions, cells are provided from the upper and lower side by medium. Culture is started with a conventional medium in a CO₂ incubator until primary adhesion of cells is replaced by adherence. Depending on applied biomaterial and cell type this process can last a few hours or even days.

Regardless of whether transparent or non transparent biomaterials are used, the degree of cell adhesion can be easily registered by fluorescence microscopy after ethanol fixation and propidium iodide incubation [29]. Such a protocol illustrates fluorescent nuclei of cells adhering to a biomaterial. For example, when MDCK cells were cultured onto four different materials such as glass, polystyrene (Thermanox®), a white or black polycarbonate filter, it was observed that each of the specimens reflects an individual growth pattern of cells ranging between perfect confluent adherence and atypical dome, blister and cluster formation respectively. This result illustrates that selected biomaterials have an enormous influence on the spatial growth pattern of adherent cells.

Not only single cells but also a thin slice of tissue can be mounted between two pieces of a woven net or a fleece within a MINUSHEET® tissue carrier [30,31]. In addition, flexible collagen sheets for adhesion of cells can be stretched like the skin of a drum in a slightly modified carrier [32,33]. The few examples show that principally a broad spectrum of biomaterials can be




Figure 1 Schematic illustration of mounting a MINUSHEET® tissue carrier. (a) First an individual biomaterial supporting optimal tissue development is selected. (b) For mounting the biomaterial measuring 13 mm in diameter is placed in the black base part of a tissue carrier. (c) After fixation by a white tension ring the mounted tissue carrier can be used for cell seeding. For easy access of medium at the basal side a tissue carrier rests with protrusions on the bottom of a dish. (d) Tissue carriers including different biomaterials are demonstrated within a 24 well culture plate.

inserted in a tissue carrier so that cells obtain an adequate surface for adhesion or adherence to develop further into a specialized tissue.

Compatible perfusion culture containers

It is obvious that the static environment within a culture dish leads to an uncontrollable increase of metabolites and a decrease of nutrition during time. Due to this reason a MINUSHEET® tissue carrier with adherent cells is used only for the relatively short period of cell seeding in the static environment of a 24 well culture plate. To provide constant nutrition and respiratory gas during further tissue development, the tissue carrier is consequently inserted in a perfusion culture container. The exact geometrical placement of the tissue carrier within a perfusion culture container guarantees on the one hand an equal distribution and on the other hand a constant provision with always fresh culture medium. The basic version of a container allows continuous medium transport to developing tissues from all sides (Figure 2a). In a gradient container the tissue carrier is held between the base and the lid so that the luminal and basal sides can be provided with individual media mimicking a typical environment for epithelia (Figure 2b). A further perfusion container is made of a transparent lid and base allowing the microscopic observation during tissue development (Figure 2c). Yet another specific container exhibits a flexible silicone lid. Applying force to this lid by an eccentric rotor simulates a mechanical load, which is required in cartilage and bone tissue engineering. Shaped tissues such as an auricle or different forms of cartilage implants can be generated with individual scaffolds in a special tissue engineering container. Finally, spatial extension of tubules derived from renal stem/progenitor cells is obtained within a perfusion culture container filled with an artificial interstitium consisting of a polyester fleece [30,31].

In conclusion, compared to other perfusion culture experiments the presently described system exhibits up to date unique features for the generation of specialized tissues. Damage during culture and transfer of growing tissue is prevented by applied MINUSHEET® tissue carriers (Figure 1). Adhesion of cells to natural extracellular matrix or synthetic biomaterials can be modulated by mounting individual supports (Figure 1a-c). While seeding of cells and initial culture can be performed in the static environment of a 24 well culture plate (Figure 1d), the exact positioning of a carrier within a perfusion container makes it possible to provide developing tissue with always fresh medium (Figure 2a). It prevents an overshoot of paracrine factors and keeps concentration of harmful metabolites low. Finally, offering different media at the luminal and basal side within a gradient culture container allows simulation of a natural barrier




Figure 2 Application of a Minusheet® tissue carrier in compatible perfusion culture containers. (a) A first type of perfusion culture container can hold six tissue carriers to provide them with always fresh medium. (b) In a gradient perfusion culture container the respective tissue is exposed to different fluids at the luminal and basal side. (c) In a microscope container with a transparent lid and base developing tissue can be observed. (d) Perfusion culture set-up is running on a laboratory table. A thermo plate maintains the desired temperature of 37°C. During culture a peristaltic pump transports the medium (1.25 ml/h) from a storage bottle (left side) to the waste bottle (right side).

environment for epithelia (Figure 2b). Thus, the described technique bridges consequently static culture in a conventional 24 well plate (Figure 1d) with modern perfusion culture in different perfusion culture containers (Figure 2).

Permanent supply of fresh culture medium

To maintain the desired temperature of 37°C within a perfusion culture container a heating plate (MEDAX-Nagel, Kiel, Germany) and a cover lid (not shown) is used (Figure 2d). The transport of culture medium is best accomplished using a slowly rotating peristaltic pump (ISMATEC, IPC N8, Wertheim, Germany), which is able to provide adjustable pump rates of 0.1 to 5 ml per hour. In personal experiments medium was transported with 1.25 ml/h for a period of at least 13 days. Adherent cells growing on a mounted tissue carrier are supplied with always fresh medium flowing from the storage bottle to the container. Thus, during long term culture the growing tissue is exposed to always fresh medium preventing in parallel an un-physiological accumulation of metabolic products and an overshoot of paracrine factors. To maintain this controlled environment for the whole culture period, the metabolized medium is not re-circulated but collected in a separate waste bottle.

Keeping the pH constant in culture medium

Working with perfusion culture equipment under atmospheric air on a laboratory table facilitates the complete handling. However, in the majority of cases conventional cell and tissue cultures are performed in a CO₂ incubator. In this case the applied media are buffered by a system consisting of a relatively high amount of NaHCO₃ and a 5% CO₂ atmosphere to maintain a constant pH between 7.2 and 7.4. If such a formulated medium is used for perfusion culture outside a CO₂ incubator, the pH will shift from the physiological range to much more alkaline values due to the low content of CO₂ (0,3%) in atmospheric air.

For that reason any medium used for perfusion culture outside a CO₂ incubator has to be stabilized by reducing the NaHCO₃ concentration and/or by adding biological buffers such as HEPES (GIBCO/Invitrogen, Karlsruhe, Germany) or BUFFER ALL (Sigma-Aldrich-Chemie, München, Germany). The necessary amount can be determined by admixing increasing amounts of biological buffer solution to an aliquot of medium. Then the medium must equilibrate over night on a thermo plate at 37°C under atmospheric air. For example, application of 50 mmol/l HEPES or an equivalent of BUFFER ALL (ca. 1%) to IMDM (Iscove's Modified Dulbecco's Medium, GIBCO/Invitrogen) maintains a constant pH of 7.4 throughout long term perfusion culture under atmospheric air on a laboratory table.

Oxygen content in transported culture medium

When a perfusion culture experiment is performed, the pH must be adjusted and sufficient respiratory gas must be present. To enrich oxygen (O_2) for example in IMDM, the medium has to be transported through a gas-permeable silicone tube. Such a silicone tube provides a large surface for gas exchange by diffusion due to a thin wall (1 mm), the small inner diameter (1 mm) and its extended length (1 m). Analysis of IMDM (3024 mg/l $NaHCO_3$, 50 mmol/l HEPES) equilibrated against atmospheric air during a standard perfusion culture experiment consequently showed partial pressures of 160 mmHg O_2 and 12 mmHg CO_2 [34].

Furthermore it has to be considered that depending on the tissue types the requirements of oxygen are individual. Therefore it is important to note that the content of oxygen can be adapted in individual perfusion culture set-ups. The technical solution is a gas exchanger module containing a gas inlet and outlet. Moreover a spiral with a long thin-walled silicone tube for medium transport is mounted inside the module. Since the tube is highly gas-permeable, optimal diffusion of gases between culture medium and surrounding atmosphere within the gas exchange module is guaranteed. In consequence, the desired gas atmosphere can be adjusted by a constant flow of a specific gas mixture through the module. That way the content of oxygen or any other gases can be modulated in the medium by diffusion. By applying this simple protocol it became possible to decrease the oxygen partial pressure within the transported medium during long term culture under absolutely sterile conditions [34]. These kinds of experiments elicited that a decrease in oxygen concentration was leading to a massive alteration of functional protein expression in generated renal collecting duct epithelia.

Elimination of harmful gas bubbles

During perfusion culture experiments gas bubbles can arise and influence in turn the flow of medium. Formation of gas bubbles is not only observed during suction of medium from the storage bottle but also during transport at material transitions between tubes and connectors.

First, these gas bubbles are so small that they cannot be recognized, but during ongoing transport of culture medium they increase in size and are able to form an embolus that massively impedes medium flow. Second, gas bubbles can accumulate in the culture container leading to a regional shortage of medium supply and breaks in the fluid continuum so that remarkable fluid pressure changes result. Most importantly, in a gradient perfusion culture experiment, in which two media are transported at exactly the same speed embolic effects can lead to pressure differences, which in turn destroy the barrier function of the contained epithelia [35].

To minimize the content of gas bubbles within a perfusion culture set-up, a gas expander module was developed [36]. This module removes gas bubbles from the medium. When medium is entering the module, it rises within a small reservoir and expands before it drops down onto a barrier. During this process gas bubbles are separated from the medium and collected at the top of the gas expander module. As a result, medium leaving the gas expander module is oxygen-saturated but free of gas bubbles.

Wide spectrum of applications

When tissues are generated in combination with new biomaterials, it is expected that they develop a high degree of functionality similar to the one known from the respective organism. Such a tissue-specific quality cannot be obtained by the use of a conventional dish, but is reached by offering an optimal surface for cell adhesion, compatible perfusion culture containers, permanent transport of nutrition including respiratory gas and constant elimination of metabolic products. Years ago such a compatible system was not available on the market. In consequence, necessary tools such as the previously described tissue carriers and versatile perfusion culture containers were developed for performance of culture experiments in the own laboratory (Figure 1 and 2).

After publication of successful results also other groups were interested in application of the innovative technique. As a result numerous papers were published dealing with the MINUSHEET® perfusion culture technique. A list of those different culture set-ups is found in the data bank 'Proceedings in perfusion culture': <http://www.biologie.uni-regensburg.de/Anatomie/Minuth/proceedings.htm>

Creating an improved fluid environment for epithelia

Experiments have shown that epithelial cells often do not develop expected cell biological features when they are cultured at the bottom of a conventional dish. Consequently, MINUSHEET® tissue carriers including an optimal biomaterial and compatible perfusion culture containers were applied to offer an improved environment.

Keeping specialized epithelia in perfusion culture

Pilot experiments were made with collecting duct (CD) cells derived from the embryonic parenchyma of neonatal kidney [37,38]. By applying the mentioned tissue carriers and perfusion culture containers polarized epithelia were harvested for the first time expressing under *in vitro* conditions cell biological features such as observed in adult Principal (P) and Intercalated Cells (IC) of the renal collecting duct [33,39-42]. In the following experiments special attention was directed to the selection of individual biomaterials improving adhesion

of cells so that they can stand fluid flow in perfusion culture [29]. Further experiments with cells derived from embryonic parenchyma brought new insights into the spatial development of renal microvasculature and glomeruli [43-46]. In addition, perfusion culture was performed to accelerate regeneration by engineered microvessels [47,48] and to investigate restoration of endothelium [49,50].

To bridge time before an implantation is made, perspectives of living conservation of human gingival epithelium was investigated in long term perfusion culture experiments [51,52]. Factors affecting reproductive aging and the development of fertilized eggs were investigated with anterior pituitary gland [53], oviduct epithelium [54] and endometrium [55]. Optimal matrix coating and adaption of continuous medium flow were elaborated for hepatocytes [56-59]. Regeneration of urothelium was investigated in combination with newly developed stent materials [60]. Effects of newly developed pharmaceuticals on ciliary beat frequency (CBF) were elaborated by differentiated nasal epithelium generated in described perfusion culture [61].

Exposure of epithelia to a luminal-basal fluid gradient

When performing perfusion culture experiments one has to consider that the environment of epithelia is changing during development. In the fetal period epithelia are exposed to the same fluid composition at the luminal and basal sides due to the still leaky barrier. However, in maturing epithelia a tight junction complex and up-regulated transport features are leading to a functional barrier. Because of the physiological seal different media are found at the luminal and basal sides [15].

When epithelial cells are cultured at the bottom of a conventional dish, all sides are exposed to the same medium. This untypical situation is leading to a permanent biological short circuit promoting cell proliferation but suppressing polar differentiation [62-67]. To overcome this problem, MINUSHEET® tissue carriers can be mounted in a gradient perfusion culture container to mimic at the luminal and basal sides a tissue-specific environment for epithelia. Following this strategy finally functional renal membranes were generated [68,69].

Originally experiments with a gradient perfusion culture container were performed to investigate the influence of different fluids at the luminal and basal sides of embryonic renal collecting duct (CD) epithelia [32,33,38,70-76]. These experiments showed that the development of a CD epithelium has an unexpected long latent period of three days. Moreover it takes at least 10 days until typical features of differentiation are up-regulated. The development is also depending on increasing NaCl concentrations specifically offered at the luminal

side (Figure 3). Surprisingly, in a continuous electrolyte gradient over many days typical epithelial features such as binding of lectins or expression of site-specific proteins such as Na/K ATPase α 5 (Figure 3a), TROMA I (Cytokeratin Endo-A; Figure 3b), cingulin (Figure 3c) and COX 2 (Figure 3d) were found to be up-regulated. However, when the increased NaCl content at the luminal side was replaced against a low NaCl concentration, gained features were found to be down-regulated. This result illustrates that the luminal-basal electrolyte gradient maintains functional features within renal epithelia.

Further challenging experiments were performed with hydrogel in MINUSHEET® tissue carriers mimicking the glomerular basement membrane. In this experimental set-up it was possible to seed endothelial cells on the one side and podocytes on the other side. Mounting these co-cultures in a gradient perfusion container, development of urine-blood barrier functions were tested




Figure 3 Cell biological features of an embryonic renal collecting duct (CD) epithelium kept for 13 days in a fluid gradient. In the gradient container at the luminal side IMDM + aldosterone (1×10^{-7} M) + 15 mmol/l NaCl, while at the basal side IMDM + aldosterone (1×10^{-7} M) was transported.

Immunohistochemistry shows a positive label for tissue-specific (a) Na/K ATPase α 5, (b) TROMA I, (c) cingulin and (d) COX 2 protein. In contrast, perfusion culture of epithelia with IMDM at the luminal and basal side does not develop these features. Basal lamina (asterisk), lumen (arrow).

[77]. In so far gradient perfusion culture made it possible to keep renal epithelia in vitro as observed under natural conditions. Furthermore such a tissue-specific development made it apparent that the applied media had to be adapted to the physiological requirements by addition of defined electrolytes [34,78,79].

Maintenance of retina

It has been shown that the retina is a complex tissue lined by a pigment epithelium that cannot be maintained in the static environment of a culture dish over prolonged periods of time. To improve the environment, intact retina was mounted onto a tissue carrier for culture within a gradient perfusion container [80-86]. It was demonstrated that pigment epithelium and neighboring neurons maintain a perfect morphology for at least 10 days. These exiting findings illustrate innovative perspectives for safety testing of newly developed pharmaceuticals designed for the intraocular medical application. Furthermore these experiments show challenging options to investigate the wide field of retina aging, degeneration and repair under realistic culture conditions [87-91].

Blood-retina and blood-brain barrier in perfusion culture

Both blood-retina and blood-brain barrier are crucial for the transport of pharmaceuticals. Gradient perfusion culture appears as an ideal technique to simulate such a blood-retina and blood-brain barrier under realistic in vitro conditions [92-94]. Indeed, the perfusion culture experiments showed new features of permeation and displayed an intact polarized expression of efflux pumps such as multidrug resistance protein (P-gp) and multi-drug resistance-associated protein (MRP).

Creating a blood-air barrier

A specific environment for pneumocytes in form of a blood-air barrier was simulated in gradient perfusion culture [28]. When pneumocytes and endothelial cells were co-cultured on a polycarbonate filter within a gradient perfusion container, development of a tight junction complex was observed sealing perfectly the blood-air barrier. Further characteristic features of polar differentiation within the epithelia were found to be up-regulated. It was argued that gradient perfusion culture in combination with pneumocytes and endothelial cells is a realistic model to investigate dose-controlled exposure of airborne particles. Moreover, to elucidate barrier transport and repair mechanisms after alveolar injury a dose controlled air-liquid interface (ALI) was investigated in a gradient perfusion container by using A549 cells [95-97].

Constructing a blood-gas barrier

Culture experiments on fish swim bladder gas gland were successfully performed in a gradient perfusion container [98]. For these series of experiments cells of gas gland were cultured on a filter at the interface between gas and culture medium. The generated epithelium revealed a typical polarity and functionality as it was observed in the environment of swim bladder gas gland in the fish.

Pharmaceutical applications applying a gradient container

Most of the administered pharmaceuticals have to pass an epithelial barrier in the organism. To test the transport of newly developed drugs through an epithelial cell layer long term gradient perfusion culture experiments were performed [99,100]. In these series of experiments it was found that Caco-2 cells are forming a tightly sealing epithelial cell layer [101,102]. Further gradient perfusion culture exhibited that reproducible results can be achieved much earlier than observed in traditional 21 day static cultures. Interestingly the permeability coefficient of several model pharmaceuticals across a Caco-2 cell layer was approximately twofold higher than obtained under conventional static conditions.

Renewal of skin

The regeneration and repair of skin is an important research area in actual biomedicine. Consequently, epidermis equivalents were generated by using gradient perfusion culture [103]. In this investigation composite grafts of INTEGRA® matrix and human keratinocytes were cultured in a gradient container in order to evaluate the potential of the cost-effective engineering of full-thickness skin grafts and the treatment of ulcers. Furthermore it was demonstrated that generation of gingival epithelium [52,104] or co-culture of keratinocytes with osteoblast-like cells kept in perfusion culture results in a much better tissue generation than observed under static culture conditions [105].

Offering an artificial interstitium for regenerating parenchyma

The repair of renal parenchyma is of special interest in acute and chronic kidney failure. To investigate the spatial regeneration of renal tubules under realistic conditions, an advanced MINUSHEET® perfusion culture protocol was developed (Figure 4). To replace coating by extracellular matrix proteins renal stem/progenitor cells were mounted between layers of a polyester fleece. The interface between the polyester fibers creates an artificial interstitium promoting spatial extension of tubules [106,107]. For example, SBA lectin labeling of a specimen kept for 13 days in perfusion culture illustrates that numerous tubules can be detected in this specific area (Figure 4a). Each of the




Figure 4 Generation of renal tubules at the interface of an artificial interstitium. (a) SBA label visualizes in an individual sample tubules developing within an artificial interstitium consisting of a polyester fleece. (b) Scanning electron microscopy illuminates growth of renal tubules (T) between fibers of the polyester fleece (PF). (c) Semithin section after Richardson staining shows generated tubules in oblique and vertical view. (d) Transmission electron microscopy demonstrates that regenerating tubules contain a polarized epithelium. Neighboring cells are separated by a tight junctional complex (arrow head). Basal lamina (asterisk), lumen (arrow).

tubules exhibits a lumen and a basal lamina. Scanning electron microscopy further illuminates the spatial growth of generated tubules between polyester fibers (Figure 4b). Semithin sections (Figure 4c) and transmission electron microscopy (Figure 4d) clearly demonstrate that generated tubules contain a polarized epithelium. In the meantime numerous data demonstrated that perfusion culture together with an artificial interstitium is an effective technique to investigate the process of tubule regeneration [108-110]. In addition, for the first time it was shown that tubule formation can be induced by aldosterone, while antagonists such as spironolactone or canrenoate prevent tubulogenic development [111-113]. By preventing the molecular contact between the mineralocorticoid receptor (MR) and heat shock protein 90 geldanamycin leads to a lack of tubule formation. In consequence, the results demonstrate that the signal initiated by aldosterone is not unspecific but is mediated via MR and depends on intact molecular interactions with heat shock proteins.

It is a new aspect that the interface between layers of polyester fleece promotes the spatial development of numerous tubules [114,115]. When such an artificial interstitium is applied, the surrounding of generating tubules is not stacked by coated extracellular matrix proteins. In consequence, it became possible for the first time to investigate the spatial growth of tubules by scanning electron microscopy and the synthesis of special interstitial molecules during regeneration (Figure 4b). Recent data further illustrate that developing tubules avoid a direct contact to each other by keeping a discrete distance during spatial development. The separation seems to be caused by an interaction between the basal lamina of generating tubules, newly synthesized fibers of the extracellular matrix and neighboring polyester fibers of the fleece [116-122].

Finally, perfusion culture experiments exhibited that application of different kinds of polyester fleeces are leading to various patterns of spatial tubule development resulting in new challenging perspectives for the regeneration of diseased renal parenchyma [31]. By performing these pilot experiments it became evident that aldosterone has a stimulating influence on the formation of renal tubules, while glucocorticoids are inducing atypical structures in form of extended cell clusters [122]. In conclusion, it is an important new finding that steroid hormones occurring in the interstitial fluid may exhibit a harmful influence on the regeneration of renal parenchyma.

Engineering of connective tissue in perfusion culture

A big research field for tissue engineering and biomaterial research is the generation of connective tissue constructs. In this coherence it is barely considered that

beside the epithelial cell layers also connective tissue can build up important barrier functions.

Connective tissue barrier

For example, experiments related to such a non-epithelial barrier were performed with dentin discs mounted in a MINUSHEET® gradient perfusion container [123-132]. Most interestingly, these data reveal that polymerized dental resin materials release residual monomers that may interact with pulp tissue. In consequence, to obtain information about diffusion of molecules through a dentin disc gradient perfusion culture appears to be an ideal model to investigate long term toxic effects under realistic in vitro conditions. Finally, new aspects of testing permeability and degradation in gelatine membranes were obtained by keeping fibroblasts in a gradient perfusion container [133].

Generation of hyaline cartilage

A great challenge in tissue engineering is the treatment of cartilage defects by regenerating chondrocytes growing onto innovative scaffold materials. In numerous cases it was demonstrated that MINUSHEET® perfusion culture can improve the cell biological quality of generated cartilage. An enormous advantage of the described technique is that damage of developing tissue within the scaffold is minimized due to the permanent elimination of biodegraded molecules by transported culture medium. For example, perfusion culture was applied successfully for the controlled regeneration of hyaline cartilage [18,21,134,135]. In these experiments it became possible to elaborate exact data concerning kinetics of the degradation process from different scaffold materials [136,137]. Furthermore the cell biological quality of generated cartilage was improved by stepwise chemical modifications of the scaffold material by perfusion culture. Following this strategy the risk of tissue repulsion after implantation could be decreased by selecting optimal scaffold materials [138-145]. Surprisingly, it was demonstrated that the application of natural extracellular matrix such as a collagen sponge does not improve the quality of generated cartilage [146]. In contrast, scaffold materials with modified polyethylene coating [147] or a gelatine-based Spongostan® [23] revealed much more cartilage specific features than observed without surface treatment. It was further detected that synovial fibroblasts are able to modulate articular matrix synthesis [148]. Finally, tissue engineering of cartilage constructs by perfusion culture seemed to be an ideal model to investigate parameters affecting destructive joint diseases [149-151].

Formation of bone constructs

Beside cartilage formation MINUSHEET® perfusion culture technique was also applied for bone tissue engineering to investigate the development of osteoblasts on ceramic

materials [24-26,152,153], decellularized spongeous bone [154], collagen membranes [155], mineralized collagen [156,157], hydroxyapatite scaffolds [158-161], PLGA sheets [162], iron based metals [163], bioactive glass [164], textile chitosan [165-167] or 3D biphasic calcium phosphate [168] scaffolds and biocorrodible bone replacement materials [169]. Finally, most important for clinical applications are experiments, which exhibit that bone development can be influenced by the sterilization procedure of scaffolds consisting of poly-d,l-lactic-co-glycolic acid [162].

A permanently occurring problem in bone tissue engineering is that unstirred and consequently harmful layers of fluid within growing tissue can arise. To compensate for this the continuous provision with nutrition and oxygen must be replaced by transportation of fluid in pulses so that generation of bone constructs results in an increased cell biological quality [170]. Last but not the least, the insights gained by bone formation may lead to an effective strategy for the regeneration of dentin [127].

Growth of muscular tissue

Surprisingly, only three papers were found that deal with the regeneration of muscular tissue in MINUSHEET® perfusion culture. When gastric mucosa was kept in a culture container, it was recognized that also smooth muscular tissue is developing within the lamina propria [67]. By applying improved biomaterials for cell seeding numerous cerebral pericytes were found to express site-specific pericytic aminopeptidase N/pAPN [171]. Proliferation and adhesion of smooth muscle cells was investigated on electrospun polymer scaffolds [172].

Development of nervous tissue

One of the main subjects in the area of neurology is the escape of dopamine synthesis during Parkinson's disease. To investigate influences affecting synthesis of dopamine, MINUSHEET® perfusion culture was performed successfully using mesencephalic neurons [173]. These experiments demonstrated for example that neurotrophins stimulate the release of dopamine via Trk and p75Lnr receptors. Furthermore it was shown in hippocampal neurons and the pheochromacytoma cell line PC 12 that application of exogenous neurotrophins exhibits positive feedback effects on secretion of synthesized neurotrophins. This pathway is mediated via an activation of tyrosine kinase neurotrophin receptors [174]. The influence of sodium in an activity-dependent secretion of neurotrophins plays thereby an important role [175]. Further significant differences in the secretion between nerve growth factor and brain-derived neurotrophic factor were observed [176]. Finally, SH-SY5Y human neuroblastoma cells were found to differentiate into a neuronal-like state in long term

perfusion culture [177]. The cells could be maintained in an active state for more than two months without the need of passaging them. In a further set of experiments RAT-1 fibroblasts were investigated expressing Cypri-dina noctiluca luciferase (CLuc) driven by the promoter of the circadian clock gene Mma11 [178]. The performed experiments showed that the CLuc reporter assay in combination with the described perfusion culture appears to be an innovative pharmacological tool for drug discovery. Finally, most promising results were obtained with fish pituitary explants to investigate vaso-tocin and isotocin release when kept in perfusion culture [179].

Conclusions

MINUSHEET® perfusion culture technique was developed to improve the environment of adherent cells and developing tissues. For optimal adhesion of cells an individual biomaterial is selected and mounted onto a tissue carrier. Seeding of cells is performed in the static environment of a conventional 24 well culture plate, while perfusion culture is performed with compatible containers offering an adequate fluid environment for the generation of specialized tissues. Numerous papers illustrate that the modular system is generating cells and tissues in excellent cell biological quality urgently needed in tissue engineering, biomaterial research and advanced pharmaceutical drug testing.

Final remarks

Numerous patents (not listed) demonstrate that Will W. Minuth has invented the presented MINUSHEET® perfusion culture technique. To introduce the developed tools on the market Katharina Lorenz-Minuth founded non-profit orientated Minucells and Minutissue GmbH (D-93077 Bad Abbach/Germany, www.minucells.com). In 1992 the project received the Philip Morris research award 'Challenge of the Future' in München/Germany.

Competing interests

The authors declare that they have no competing interests.

Authors' contributions

WM contributed to the collection of information, interpretation of data and writing of the manuscript. LD performed the generation of renal tubules, histochemistry, electron microscopy, design of figures and critical revision of the manuscript. All authors read and approved the final manuscript.

Received: 27 February 2012 Accepted: 2 October 2012

Published: 5 October 2012

References

- Elaut G, Henkens T, Papeleu P, Snykers S, Vinken M, Vanhaecke T, Rogiers V: Molecular mechanisms underlying the dedifferentiation process of isolated hepatocytes and their cultures. *Curr Drug Metab* 2006, **7**:629–660.
- Schuh E, Hofmann S, Stok K, Notbohm H, Müller R, Rotter N: Chondrocyte redifferentiation in 3D: the effect of adhesion site density and substrate elasticity. *J Biomed Mater Res A* 2011, **100**:38–47.
- Zhang Y, Li TS, Lee ST, Wawrowsky KA, Cheng K, Galang G, Malliaras K, Abraham MR, Wang C, Marban E: Dedifferentiation and proliferation of mammalian cardiomyocytes. *PLoS One* 2010, **5**e12559.
- Liu Y, Jiang X, Yu MK, Dong J, Zhang X, Tsang LL, Chung YW, Li T, Chan HC: Switching from bone marrow-derived neurons to epithelial cells through dedifferentiation and translineage redifferentiation. *Cell Biol Int* 2010, **34**:1075–1083.
- Russ HA, Bar Y, Ravassard P, Efrat S: In vitro proliferation of cells derived from adult human beta-cells revealed by cell-lineage tracing. *Diabetes* 2008, **6**:1575–1583.
- Flieger A, Golka K, Schulze H, Follmann W: Primary cultures of human urothelial cells for genotoxicity testing. *J Toxicol Environ Health A* 2008, **71**:930–935.
- Cournil-Henrionnet C, Huselstein C, Wang Y, Galois L, Mainard D, Decot V, et al: Phenotypic analysis of cell surface markers and gene expression of human mesenchymal stem cells and chondrocytes during monolayer expansion. *Biorheology* 2008, **45**:513–526.
- Stoddart MJ, Grad S, Eglin D, Alini M: Cells and biomaterials in cartilage tissue engineering. *Regen Med* 2009, **4**:81–98.
- Mikos AG, Herring SW, Ochareon P, Elisseeff J, Lu HH, Kandel R, et al: Engineering complex tissues. *Tissue Eng* 2006, **12**:3307–3339.
- Sands RW, Mooney DJ: Polymers to direct cell fate by controlling the microenvironment. *Curr Opin Biotechnol* 2007, **18**:448–453.
- Hwang NS, Varghese S, Elisseeff J: Controlled differentiation of stem cells. *Adv Drug Deliv Rev* 2008, **60**:199–214.
- Dawson JL, Wahl DA, Lanham SA, Kanczler JM, Czernuszka JT, Oreffo ROC: Development of specific collagen scaffolds to support the osteogenic and chondrogenic differentiation of human bone marrow stromal cells. *Biomaterials* 2008, **29**:3105–3116.
- Lopez JL, Mouw JK, Weaver VM: Biomechanical regulation of cell orientation and fate. *Oncogene* 2008, **27**:6981–6993.
- Potter W, Kalil RE, Kao WJ: Biomimetic material systems for neural progenitor cell-based therapy. *Front Biosci* 2008, **13**:806–821.
- Minuth WW, Denk L: Advanced culture experiments with adherent cells. From single cells to specialized tissues in perfusion culture: University of Regensburg; 2011. URN:nbn:de:bvb:355-epub-214847. ISBN 978-3-88246-330-9.
- Minuth WW, Denk L, Glashausser A: A modular culture system for the generation of multiple specialized tissues. *Biomaterials* 2010, **31**:2945–2954.
- von Recum HA, Cleek RL, Eskin SG, Mikos AG: Degradation of polydispersed poly(L-lactic acid) to modulate lactic acid release. *Biomaterials* 1995, **16**:441–447.
- Thomson RC, Yaszemski MJ, Powers JM, Mikos AG: Fabrication of biodegradable polymer scaffolds to engineer trabecular bone. *J Biomater Sci-Polym E* 1995, **7**:23–38.
- Hutmacher D, Hürzeler MB, Schliephake H: A review of material properties of biodegradable and bioresorbable polymers and devices for GTR and GBR applications. *Int J Oral Maxillofac Implants* 1996, **11**:667–678.
- Bujia J, Sittoner M, Minuth WW, Hammer C, Burmester G, Kastenbauer E: Engineering of cartilage tissue using bioresorbable polymer fleeces and perfusion culture. *Acta Otolaryngol* 1995, **115**:307–310.
- Sittoner M, Bujia J, Minuth WW, Hammer C, Burmester GR: Engineering of cartilage tissue using bioresorbable polymer carriers in perfusion culture. *Biomaterials* 1994, **15**:451–456.
- Chen DC, Avansino JR, Agopian VG, Hoagland VD, Woolman JD, Pan S, et al: Comparison of polyester scaffolds for bioengineered intestinal mucosa. *Cells Tissues Organs* 2006, **184**:154–165.
- Anders JO, Mollenhauer J, Beberhold A, Kinne RW, Venbrooks RA: Gelatin-based haemostyptic Spongostan as a possible three-dimensional scaffold for a chondrocyte matrix? An experimental study with bovine chondrocytes. *J Bone Joint Surg Br* 2009, **91**:409–416.
- Uemura T, Dong J, Wang Y, Kojima H, Saito T, Iejima D, et al: Transplantation of cultured bone cells using combinations of scaffolds and culture techniques. *Biomaterials* 2003, **24**:2277–2286.
- Wang Y, Uemura T, Dong J, Kojima H, Tanaka J, Tateishi T: Application of perfusion culture system improves in vitro and in vivo osteogenesis of bone marrow-derived osteoblastic cells in porous ceramic materials. *Tissue Eng* 2003, **9**:1205–1214.

26. Leukers B, Gürkan H, Irsen SH, Milz S, Tille C, Seitz H, et al: Biocompatibility of ceramic scaffolds for bone replacement made by 3D printing. *Materialwiss Werkst* 2005, **36**:781–787.
27. Gueven N, Glatthaar B, Manke HG, Haemmerle H: Co-cultivation of rat pneumocytes and bovine endothelial cells on a liquid-air interface. *Eur Respir J* 1996, **9**:968–975.
28. Gläum R, Wiedmann-Al-Ahmad M, Huebner U, Schmelzeisen R: Tissue engineering of composite grafts: cocultivation of human oral keratinocytes and human osteoblast-like cells on laminin-coated polycarbonate membranes and equine collagen membranes under different culture conditions. *J Biomed Mater Res* 2009, **93**:704–715.
29. Minuth WW, Majer V, Kloth S, Dermietzel R: Growth of MDCK cells on non-transparent supports. *In Vitro Cell Dev Biol Anim* 1994, **30**:12–14.
30. Minuth WW, Denk L, Heber S: Growth of embryonic renal parenchyme at the interphase of a polyester artificial interstitium. *Biomaterials* 2005, **26**:6588–6598.
31. Roessger A, Denk L, Minuth WW: Potential of stem/progenitor cell cultures within polyester fleeces to regenerate renal tubules. *Biomaterials* 2009, **30**:723–732.
32. Minuth WW, Dermietzel R, Kloth S, Henerkes B: A new method culturing renal cells under permanent perfusion and producing a luminal-basal medium gradient. *Kidney Int* 1992, **41**:215–219.
33. Minuth WW, Aigner J, Kloth S, Steiner P, Tauc M, Jennings ML: Culture of embryonic renal collecting duct epithelia kept in a gradient container. *Ped Nephrology* 1997, **11**:140–147.
34. Strehl R, Schumacher K, Minuth WW: Controlled respiratory gas delivery to embryonic renal epithelial explants in perfusion culture. *Tissue Eng* 2004, **10**:1196–1203.
35. Minuth WW, Denk L, Roessger A: Gradient perfusion culture-simulating a tissue-specific environment for epithelia in biomedicine. *J Epithel Biology & Pharmacology* 2009, **2**:1–13.
36. Minuth WW, Strehl R, Schumacher K: Tissue factory: conceptual design of a modular system for the in vitro generation of functional tissues. *Tissue Eng* 2004, **10**:285–294.
37. Minuth WW, Rudolph U: A compatible support system for cell culture in biomedical research. *Cytotechnology* 1990, **4**:181–189.
38. Minuth WW, Stöckl G, Kloth S, Demietzel R: Construction of an apparatus for cell and tissue cultures which enables in vitro experiments under organo-typical conditions. *Eur J Cell Biol* 1992, **57**:132–137.
39. Herter P, Laube G, Gronczewski J, Minuth WW: Silver-enhanced colloidal-gold labelling of rabbit kidney collecting-duct cell surfaces imaged by scanning electron microscopy. *J Microsc* 1993, **171**:107–115.
40. Minuth WW, Fietzek W, Kloth S, Aigner J, Herter P, Röckl W, et al: Aldosterone modulates PNA binding cell isoforms within renal collecting duct epithelium. *Kidney Int* 1993, **44**:537–544.
41. Aigner J, Kloth S, Kubitzka M, Kashgarian M, Dermietzel R, Minuth WW: Maturation of renal collecting duct cells in vivo and under perfusion culture. *Epithelial Cell Biol* 1994, **3**:70–78.
42. Aigner J, Kloth S, Jennings ML, Minuth WW: Transitional differentiation patterns of principal and intercalated cells during renal collecting duct development. *Epithelial Cell Biol* 1995, **4**:121–130.
43. Kloth S, Schmidbauer A, Kubitzka M, Weich HA, Minuth WW: Developing renal microvasculature can be maintained under perfusion culture conditions. *Eur J Cell Biol* 1994, **63**:84–95.
44. Kloth S, Ebenbeck C, Kubitzka M, Schmidbauer A, Röckl W, Minuth WW: Stimulation of renal microvascular development under organotypic culture conditions. *FASEB J* 1995, **9**:963–967.
45. Kloth S, Gerdes J, Wanke C, Minuth WW: Basic fibroblast growth factor is a morphogenic modulator in kidney vessel development. *Kidney Int* 1998, **53**:970–978.
46. Kloth S, Suter-Cazzolara C: Modulation of renal blood vessel formation by glial cell line-derived neurotrophic factor. *Microvasc Res* 2000, **1**:190–194.
47. Frerich B, Zückmantel K, Hemprich A: Microvascular engineering in perfusion culture: immunohistochemistry and CLMS findings. *Head Face Med* 2006, **2**:26.
48. Frerich B, Zückmantel K, Winter K, Müller-Dürwald S, Hemprich A: Maturation of capillary-like structures in a tube-like construct in perfusion and rotation culture. *Int J Oral Maxillofac Surg* 2008, **37**:459–466.
49. Bakowsky U, Ehrhardt C, Loebach C, Li P, Kneuer C, Jahn D, et al: Adhesion molecule-modified cardiovascular prothesis: Characterization of cellular adhesion in a cell culture model and by cellular force spectroscopy. In *Current Research and Applications*. Weinheim: Wiley-VCH; 2005:157–173.
50. Hayashi M, Matsuzaki Y, Shimonaka M: Impact of plasminogen on an in vitro wound healing model based on a perfusion cell culture system. *Mol Cell Biochem* 2009, **322**:1–13.
51. Lehmann P, Kloth S, Aigner J, Dammer R, Minuth WW: Lebende Langzeitkonserverung von humarer Gingiva in der Perfusionskultur. *Mund Kiefer Gesichtschir* 1997, **1**:26–30.
52. Lauer G: Optimizing the mucosa graft: Developing gingival keratinocyte-fibroblast contact. In *Fundamentals of tissue engineering and regenerative medicine 2009, Kap 27.6*. Edited by Meyer U, et al. Heidelberg: Springer Verlag:375–380.
53. Zheng W, Jimenez-Linan M, Rubin BS, Halvorson LM: Anterior pituitary gene expression with reproductive aging in the female rat. *Biol Reprod* 2007, **76**:1091–1102.
54. Reischl J, Prelle K, Schöl H, Neumüller C, Einspanier R, Sinowitz F, et al: Factors affecting proliferation and dedifferentiation of primary bovine oviduct epithelial cells in vitro. *Cell Tissue Res* 1999, **296**:371–383.
55. Tiemann U, Bücher K, Pfarrer C, Pöhland R, Becker F, Kanitz W, et al: Influence of ovarian steroid hormones or platelet-activating factor on mRNA of platelet-activating factor receptor in endometrial explant perfusion cultures from ovariectomized bovine. *Prostaglandins Other Lipid Mediat* 2005, **76**:35–47.
56. Fiegel HC, Havers J, Kneser U, Smith MK, Moeller T, Kluth D, et al: Influence of flow conditions and matrix coatings on growth and differentiation of three-dimensionally cultured rat hepatocytes. *Tissue Eng* 2004, **10**:165–174.
57. Schumacher K, Khong Y, Chang S, Ni J, Sun W, Yu H: Perfusion culture improves the maintenance of cultured liver tissue slices. *Tissue Eng* 2007, **13**:197–205.
58. Du Y, Han R, Wen F, Ng San San S, Xia L, Wohland T, et al: Synthetic sandwich culture of 3D hepatocyte monolayer. *Biomaterials* 2008, **29**:290–301.
59. Xia L, Ng S, Han R, Tu X, Xiao G, Leo HL, et al: Laminar-flow immediate-overlay hepatocyte sandwich perfusion system for drug hepatotoxicity testing. *Biomaterials* 2009, **30**:5927–5936.
60. Sternberg K, Selent C, Hakanßen N, Töllner J, Langer T, Seiter H, et al: Bioartifizielle Materialien in der Urologie. *Urologe* 2004, **10**:1200–1207.
61. Dimova S, Vlaeminck V, Brewster ME, Noppe M, Jorissen M, Augustijns P: Stable ciliary activity in human nasal epithelial cells grown in a perfusion system. *Int J Pharm* 2005, **292**:157–168.
62. Bianco P, Robey PG: Stem cells in tissue engineering. *Nature* 2001, **414**:118–121.
63. Bhadriraju K, Hansen LK: Hepatocyte adhesion, growth and differentiated function on RGD-containing proteins. *Biomaterials* 2000, **21**:267–272.
64. Pilot F, Lecuit T: Compartmentalized morphogenesis in epithelia: from cell to tissue shape. *Dev Dyn* 2005, **232**:685–694.
65. Minuth WW, Aigner J, Kubat B, Kloth S: Improved differentiation of renal tubular epithelium in vitro Potential for tissue engineering. *Expt Nephrol* 1997, **5**:10–17.
66. Klug MG, Soonpaa MH, Koh GY, Field LJ: Genetically selected cardiomyocytes from differentiating embryonic stem cells form stable intracardiac grafts. *J Clin Invest* 1996, **98**:216–224.
67. Kloth S, Eckert E, Klein SJ, Monzer J, Wanke C, Minuth WW: Gastric epithelium under organotypic perfusion culture. *In Vitro Cell Dev Biol Anim* 1998, **7**:515–517.
68. Dankers PY, Boomker JM, Huizinga-van der Vlag A, Wisse E, Smets FM, Harmsen MC, Bosman AW, Meijer W, van Lyn MJ: Bioengineering of living renal membranes consisting of hierarchical, bioactive supramolecular meshes and human tubular cells. *Biomaterials* 2011, **32**:723–733.
69. Dankers PY, Boomker JM, Huizinga-van der Vlag A, Smedts FM, Harmsen MC, van Lyn MJ: The use of fibrous, supramolecular membranes and human tubular cells for renal epithelial tissue engineering towards a suitable membrane for a bioartificial kidney. *Macromol Biosci* 2010, **10**:1345–1354.
70. Minuth WW: Methode zur Kultivierung von Zellen. Patent 1990, DE-PS 3923279.
71. Minuth WW, Steiner P, Strehl R, Kloth S: Electrolyte environment modulated differentiation in embryonic renal collecting duct epithelium. *Exp Nephrol* 1997, **5**:414–422.

72. Minuth WW, Steiner P, Strehl R, Schumacher K, de Vries U, Kloth S: Modulation of cell differentiation in perfusion culture. *Exptl Nephrology* 1999, **7**:394–406.
73. Minuth WW, Strehl R, Schumacher K, de Vries U: Long term culture of epithelia in a continuous fluid gradient for biomaterial testing and tissue engineering. *J Biomat Sci Polymer Edn* 2001, **12**:353–365.
74. Minuth WW, Schumacher K, Strehl R: Renal epithelia in long term gradient culture for biomaterial testing and tissue engineering. *Biomed Mater Eng* 2005, **15**:51–63.
75. Steiner P, Strehl R, Kloth S, Tauc M, Minuth WW: In vitro development and preservation of specific features of collecting duct epithelial cells from embryonic rabbit kidney are regulated by the electrolyte environment. *Differentiation* 1997, **62**:193–202.
76. Schumacher K, Strehl R, Vries U, Minuth WW: Advanced technique for long term culture of epithelia in a continuous luminal - basal medium gradient. *Biomaterials* 2002, **23**:805–815.
77. Bruggeman LA, Doan RP, Loftis J, Darr A, Calabro A: A cell culture system for the structure and hydrogel properties of basement membranes: applications to capillary walls. *Cell Mol Bioeng* 2012, **5**:194–204.
78. Schumacher K, Strehl R, Kloth S, Tauc M, Minuth WW: The influence of culture media on embryonic renal collecting duct cell differentiation. *In Vitro Cell Dev Biol Anim* 1999, **8**:465–471.
79. Schumacher K, Castrop H, Strehl R, de Vries U, Minuth WW: Cyclooxygenases in the collecting duct of neonatal rabbit kidney. *Cell Physiol Biochem* 2002, **12**:63–74.
80. Framme C, Kobuch K, Eckert E, Monzer J, Roider J: RPE in perfusion tissue culture and its response to laser application. *Ophthalmologica* 2002, **216**:320–328.
81. Spiegel D, Scheftaler M, Kobuch K: Outflow facilities through Descemet's membrane in rabbits. *Graefes Arch Clin Exp Ophthalmol* 2002, **240**:111–113.
82. Saikia P, Maisch T, Kobuch K, Jackson TL, Bäumer W, Szeimies RM, et al: Safety testing of indocyanine green in an ex vivo porcine retina model. *Invest Ophthalmol Vis Sci* 2006, **47**:4998–5003.
83. Jian GE, Jingbo L: The stem cell and tissue engineering research in chinese ophthalmology. *Front Med China* 2007, **1**:6–10.
84. Hamilton RD, Foss AJ, Leach L: Establishment of a human in vitro model of the outer blood-retinal barrier. *J Anat* 2007, **211**:707–716.
85. Hammer M, Richter S, Kobuch K, Mata N, Schweitzer D: Intrinsic tissue fluorescence in an organotypic perfusion culture of the porcine ocular fundus exposed to blue light and free radicals. *Graefes Arch Clin Exp Ophthalmol* 2008, **246**:979–988.
86. Kobuch K, Hermann WA, Framme C, Sachs HG, Gabel VP: Maintenance of adult porcine retina and retinal pigment epithelium in perfusion culture: characterization of an organotypic in vitro model. *Exp Eye Res* 2008, **86**:661–668.
87. Klettner A, Roider J: Constitutive and oxidative-stress-induced expression of VEGF in the RPE are differently regulated by different mitogen-activated protein kinases. *Graefes Arch Clin Exp Ophthalmol* 2009, **247**:1487–1492.
88. Klettner A, Kruse M, Meyer T, Welch D, Kabelitz D, Roider J: Different properties of VEGF-antagonists: bevacizumab but not Ranibizumab accumulates in RPE cells. *Graefes Arch Clin Exp Ophthalmol* 2009, **247**:1601–8.
89. Miura Y, Klettner A, Noelle B, Hasselbach H, Roider J: Change of morphological and functional characteristics of retinal pigment epithelium cells during cultivation of retinal pigment epithelium-choroid perfusion tissue culture. *Ophthalmic Res* 2010, **43**:122–133.
90. Klettner AK, Doths J, Roider J: Nicotine reduces VEGF-secretion and phagocytic activity in porcine RPE. *Graefes Arch Clin Exp Ophthalmol* 2011, **250**:33–8.
91. Treumer F, Klettner A, Baltz J, Hassain AA, Brinkmann AA, Roider J, Hillenkamp J: Vectorial release of matrix metalloproteinases (MMPs) from porcine RPE-choroid explants following selective retina therapy (SRT): Towards slowing the macular aging process. *Exp Eye Res* 2012, **97**:63–72.
92. Steuer H, Jaworski A, Stoll D, Schlosshauer B: In vitro model of the outer blood-retina barrier. *Brain Res Protocols* 2004, **13**:26–36.
93. Steuer H, Jaworski A, Eleger B, Kaussmann M, Keldenich J, Schneider H, et al: Functional characterization and comparison of the outer blood-retina barrier and the blood-brain barrier. *Invest Ophthalmol Vis Sci* 2005, **46**:1047–1053.
94. Hamilton RD, Leach L: Isolation and properties of an in vitro human outer blood-retinal barrier model. *Methods Mol Biol* 2011, **686**:401–416.
95. Tippe A, Heinemann U, Roth C: Deposition of fine and ultrafine aerosol particles during exposure at the air/cell interface. *Aerosol Sci* 2002, **33**:2007–2218.
96. Bitterle E, Karg E, Schroepel A, Kreyling WG, Tippe A, Ferron GA, et al: Dose-controlled exposure of A549 epithelial cells at the air-liquid interphase to airborne ultrafine carbonaceous particles. *Chemosphere* 2006, **65**:1784–1790.
97. Maier KL, Alessandrini F, Beck-Speier I, Hofer TPJ, Diabaté S, Bitterle E, et al: Health effects of ambient particulate matter—biological mechanisms and inflammatory responses to in vitro and in vivo particle exposures. *Inhal Toxicol* 2008, **20**:319–337.
98. Prem C, Pelster B: Swimbladder gas gland cells of the European eel cultured in superfusion system. *Methods Cell Sci* 2000, **22**:125–132.
99. Kloth S, Kobuch K, Domokos J, Wanke C, Minuth WW: Interactive tissue culture systems: innovative tools for toxicity testing. *BIOforum Int* 1999, **3**:70–72.
100. Kloth S, Kobuch K, Domokos J, Wanke C, Monzer J: Polar application of test substances in an organotypic environment and under continuous flow: a new tissue-based test concept for a broad range of applications in pharmacotoxicity. *Toxicology in vitro* 2000, **14**:265–274.
101. Masungi C, Borremans C, Willems B, Mensch J, van Dijck A, Augustijns P, et al: Usefulness of a novel Caco-2 cell perfusion system. I. In vitro prediction of the absorption potential of passively diffused compounds. *J Pharm Sci* 2004, **93**:2507–2521.
102. Masungi C, Mensch J, Willems B, van Dijck A, Borremans C, Noppe M, et al: Usefulness of a novel Caco-2 cell perfusion system II. Characterization of monolayer properties and peptidase activity. *Die Pharmazie* 2009, **64**:36–42.
103. Kremer M, Lang E, Berger A: Organotypical engineering of differentiated composite-skin equivalents of human keratinocytes in a collagen-GAG matrix (INTEGRA Artificial Skin) in a perfusion culture system. *Langenbecks Arch Surg* 2001, **386**:357–363.
104. Hagedorn GM, Blank A, Weiland B, Spassov A, Lauer G: Perfusion culture promotes differentiation of oral keratocytes in vitro. *J Physiol Pharmacol* 2009, **60**:25–9.
105. Glauert R, Wiemann-Al-Ahmad M, Huebner U, Schmelzleisen R: Cocultivation of human keratinocytes and human osteoblast-like cells on laminin-coated polycarbonate membranes under different culture conditions. *J Biomed Mater Res A* 2009, **93**:704–714.
106. Minuth WW, Schumacher K: From the renal stem cell niche to functional tubule. *Med Klin (Munich)* 2003, **15**:31–35.
107. Minuth WW, Sorokin L, Schumacher K: Generation of renal tubules at the interface of an artificial interstitium. *Cell Physiol Biochem* 2004, **14**:387–394.
108. Heber S, Denk L, Hu K, Minuth WW: Modulating the development of renal tubules growing in serum-free culture medium at an artificial interstitium. *Tissue Eng* 2007, **13**:281–292.
109. Hu K, Denk L, de Vries U, Minuth WW: Chemically defined medium environment for the development of renal stem cells into tubules. *Biotechnol J* 2007, **2**:992–995.
110. Minuth WW, Denk L, Hu K: The role of polyester interstitium and aldosterone during structural development of renal tubules in serum-free medium. *Biomaterials* 2007, **30**:4418–4428.
111. Minuth WW, Denk L, Hu K, Castrop H, Gomez-Sanchez C: The tubulogenic effect of aldosterone is attributed to intact binding and intracellular response of the mineralocorticoid receptor. *CEJB* 2007, **2**:307–325.
112. Minuth WW, Denk L: Aldosterone-dependent generation of tubules derived from renal stem/progenitor cells. *Transplantationsmedizin* 2008, **20**:42–47.
113. Minuth WW, Blattmann A, Denk L, Castrop H: Mineralocorticoid receptor, heat shock proteins and immunophilins participate in the transmission of the tubulogenic signal of aldosterone. *J Epithel Biol Pharmacol* 2008, **11**:24–34.
114. Blattmann A, Denk L, Strehl R, Castrop H, Minuth WW: The formation of pores in the basal lamina of regenerated renal tubules. *Biomaterials* 2008, **29**:2749–2756.
115. Minuth WW, Denk L, Meese C, Rachel R, Roessiger A: Ultrastructural insights in the interface between generated renal tubules and a polyester interstitium. *Langmuir* 2009, **25**:4621–4627.

116. Glashäuser A, Denk L, Minuth WW: Polyester fleeces used as an artificial interstitium influence the spatial growth of regenerating tubules. *J Tissue Sci Eng* 2011, **2**:105.
117. Minuth WW, Denk L, Roessger A: Regenerating tubules for kidney repair. In *Stem Cell Engineering-Principals and Applications*. Edited by Artmann, Minger, Hescheler. Heidelberg: Springer Verlag; 2010:321–344.
118. Minuth WW, Denk L, Glashäuser A: Towards a guided regeneration of renal tubules at a polyester interstitium. *Materials* 2010, **3**:2369–2392.
119. Miess C, Glashäuser A, Denk L, de Vries U, Minuth WW: The interface between generating renal tubules and a polyester fleece in comparison to the interstitium of the developing kidney. *Ann Biomed Eng* 2010, **38**:2197–2209.
120. Minuth WW, Denk L, Glashäuser A: Cell and drug delivery therapeutics for renal tissue regeneration. *Adv Drug Deliv Rev* 2010, **62**:841–854.
121. Minuth WW, Denk L: Interstitial interfaces show marked differences in regenerating tubules, matured tubules, and the renal stem/progenitor cell niche. *J Biomed Mater Res A* 2012, **100**:1115–25.
122. Minuth WW, Denk L, Glashäuser A: Promoting and harmful effects of steroid hormones on renal stem/progenitor cells development. *J Tissue Sience Engeering* 2010, **1**:101.
123. Schmalz G, Garhammer P, Schweikl H: A commercially available cell culture device modified for dentin barrier test. *J Endod* 1996, **22**:249–252.
124. Schmalz G, Schuster U, Nuetzel K, Schweikl H: An *in vitro* pulp chamber with three-dimensional cell cultures. *J Endod* 1999, **25**:24–29.
125. Schmalz G, Schuster U, Thonemann B, Barth M, Esterbauer S: Dentin barrier test with transfected bovine-pulp derived cells. *J Endod* 2001, **27**:96–102.
126. Schmalz G, Schuster U, Koch A, Schweikl H: Cytotoxicity of low pH dentin-bonding agents in a dentin-barrier test *in vitro*. *J Endod* 2002, **28**:188–192.
127. Camps J, About L, Thonneman B, Mitsiadis TA, Schmaltz G, Franquin JC: Two- versus three-dimensional *in vitro* differentiation of human pulp cells into odontoblastic cells. *Connect Tissue Res* 2002, **43**:396–400.
128. Galler K, Hiller KA, Ettl T, Schmalz G: Selective influence of dentin thickness upon toxicity of dentin contacting materials. *J Endod* 2005, **31**:396–399.
129. Demirci M, Hiller KA, Bosi C, Galler K, Schmalz G, Schweikl H: The induction of oxidative stress, cytotoxicity and genotoxicity by dental adhesives. *Dent Mater* 2008, **24**:362–371.
130. Vajrabhaya L, Korsuwanawong S, Bosi C, Schmalz G: The cytotoxicity of self-etching primer bonding agents *in vitro*. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 2009, **3**:e86–90.
131. Ulker HE, Sengun A: Cytotoxicity evaluation of self adhesive composite resin cements by dentin barrier test on 3D pulp cells. *Eur J Dent* 2009, **3**:120–126.
132. Sengun A, Yalcin M, Ulker HE, Oztürk B, Hakki SS: Cytotoxicity evaluation of dentin bonding agents by dentin barrier test on 3-dimensional pulp cells. *Oral Surg Oral Pathol Oral Radiol Endond* 2011, **112**:e83–8.
133. Dreesmann L, Hajosch R, Ahlers M, Nuernberger JV, Schlosshauer B: Permeability testing of biomaterial membranes. *Biomed Mater* 2008, **3**:34119.
134. Bujia J, Sittinger M, Hammer C, Burmester G: Züchtung menschlichen Knorpelgewebes mit Hilfe einer Perfusionskammer. *Laryngo Rhino Otol* 1994, **73**:577–580.
135. Sittinger M, Bujia J, Rotter N, Reitzel D, Minuth WW, Burmester GR: Tissue engineering and autologous transplant formation: practical approaches with resorbable biomaterials and new cell culture techniques. *Biomaterials* 1996, **17**:237–242.
136. Sittinger M, Schultz O, Keyszer G, Minuth WW, Burmester GR: Artificial tissues in perfusion culture. *Int J Artif Organs* 1997, **20**:57–62.
137. Guarnizo AC, Rodriguez EV, Albaladejo JM, Carrasquer TC, Manzano JA, Cacha JP, et al: Autoinjerto subcutaneo de cartílago neosintetizado utilizando el polímero ETHISORB en conejos. *Acta Otorrinolaringol Esp* 2002, **53**:631–636.
138. Schneider T, Kohl B, Sauter T, Becker T, Kratz K, Schossig M, Hiebl B, Jung F, Lendlein A, Ertel W, Schulze-Tanzil G: Viability, adhesion and differentiated phenotype of articular chondrocytes on degradable polymers and electro-spun structures thereof. *Macromol Symp* 2011, **309**/310:28–39.
139. Rotter N, Aigner J, Naumann A, Planck H, Hammer C, Burmester G, et al: Cartilage reconstruction in head and neck surgery: comparison of resorbable polymer scaffolds for tissue engineering of human septal cartilage. *J Biomed Mater Res* 1998, **42**:347–356.
140. Kreklau B, Sittinger M, Mensing MB, Voigt C, Berger G, Burmester GR, et al: Tissue engineering of biphasic joint cartilage transplants. *Biomaterials* 1999, **20**:1743–1749.
141. Rotter N, Aigner J, Naumann A, Hammer C, Sittinger M: Behavior of tissue-engineered human cartilage after transplantation into nude mice. *J Mater Sci Mater Med* 1999, **10**:689–693.
142. Duda GN, Haisch A, Endres M, Gebert C, Schroeder D, Hoffmann JE, et al: Mechanical quality of tissue engineered cartilage: results after 6 and 12 weeks *in vivo*. *J Biomed Mater Res* 2000, **53**:673–677.
143. Haisch A, Kläring S, Gröger A, Gebert C, Sittinger M: A tissue-engineering model for the manufacture of auricular-shaped cartilage implants. *Eur Arch Otorhinolaryngol* 2002, **259**:316–321.
144. Duda GN, Kliche A, Kleemann R, Hoffmann JE, Sittinger M, Haisch A: Does low-intensity pulsed ultrasound stimulate maturation of tissue-engineered cartilage? *J Biomed Mater Res B Appl Biomater* 2004, **68**:21–28.
145. Gille J, Meisner U, Ehlers EM, Müller A, Russlies M, Behrens P: Migration pattern, morphology and viability of cells suspended in or sealed with fibrin glue: a histomorphologic study. *Tissue Cell* 2005, **37**:339–348.
146. Fuss M, Ehlers EM, Russlies M, Rohwedel J, Behrens P: Characteristics of human chondrocytes, osteoblasts and fibroblasts seeded onto a type I/III collagen sponge under different culture conditions. A light, scanning and transmission electron microscopy study. *Ann Anat* 2000, **182**:303–310.
147. Röpke E, Schön I, Vogel J, Jamali J, Bloching M, Berghaus A: Screening von modifizierten Polyethylenoberflächen für das Tissue engineering von Chondrozyten. *Laryngo-Rhino-Otol* 2007, **86**:37–43.
148. Steinhagen J, Niggemeyer O, Fuerst M, Rüther W, Schünke M, Kurz B: Perfusion culture system: Synovial fibroblasts modulare articular chondrocyte matrix synthesis *in vitro*. *Tissue Cell* 2010, **42**:151–157.
149. Schultz O, Keyszer G, Zacher J, Sittinger M, Burmester GR: Development of *in vitro* model systems for destructive joint diseases: novel strategies for establishing inflammatory pannus. *Arthritis Rheum* 1997, **40**:1420–1428.
150. Risbud MV, Sittinger M: Tissue engineering: advances in *in vitro* cartilage generation. *Trends Biotechnol* 2002, **20**:351–356.
151. Bücheler M, Haisch A: Tissue engineering in otorhinolaryngology. *DNA Cell Biol* 2003, **22**:549–564.
152. Yeatts A, Fisher JP: Bone tissue engineering bioreactors: dynamic culture and the influence of shear stress. *Bone* 2010, **48**:171–81.
153. Bernhardt A, Lode A, Peters F, Gelinsky M: Optimization of culture conditions for osteogenically-induced mesenchymal stem cells in β-tricalcium phosphate ceramics with large interconnected channels. *J Tissue Regen Med* 2010, **5**:444–53.
154. Seitz S, Ern K, Lamper G, Docheva D, Drosse I, Milz S, et al: Influence of *in vitro* cultivation on the integration of cell-matrix constructs after subcutaneous implantation. *Tissue Eng* 2007, **13**:1059–1067.
155. Rothamel D, Schwarz F, Sculean A, Herten M, Scherbaum W, Becker J: Biocompatibility of various collagen membranes in cultures of human PDL fibroblasts and human osteoblast-like cells. *Clin Oral Implants Res* 2004, **15**:443–449.
156. Gelinsky M, König U, Sewing A, Pompe W: Poröse Scaffolds aus mineralisiertem Kollagen - ein biomimetisches Knochenersatzmaterial. *Materialwiss Werkst* 2004, **35**:229–233.
157. Bernhardt A, Lode A, Boxberger S, Pompe W, Gelinsky M: Mineralised collagen—an artificial, extracellular bone matrix-improves osteogenic differentiation of bone marrow stromal cells. *J Mater Sci Mater Med* 2008, **19**:269–275.
158. Leukers B, Gülkhan H, Irseñ SH, Milz S, Tille C, Schieker M, et al: Hydroxyapatite scaffolds for bone tissue engineering made by 3D printing. *J Mater Sci Mater Med* 2005, **16**:1121–1124.
159. Detsch R, Uhl F, Deisinger U, Ziegler G: 3D-Cultivation of bone marrow stromal cells on hydroxyapatite scaffolds fabricated by dispense-plotting and negative mould technique. *J Mater Sci Mater Med* 2008, **19**:1491–1496.
160. da Silva HM, Mateescu M, Damia C, Champion E, Soares G, Anselme K: Importance of dynamic culture for evaluating osteoblast activity on dense silicon-substituted hydroxyapatite. *Colloids Surf B Biointerfaces* 2010, **80**:138–44.
161. da Silva HM, Mateescu M, Ponche A, Damia C, Champion E, Soares G, Anselme K: Surface transformation of silicon-doped hydroxyapatite immersed in culture medium under dynamic and static conditions. *Colloids Surf B Biointerfaces* 2010, **75**:349–55.

162. Shearer H, Ellis MJ, Perera SP, Chaudhuri JB: Effects of common sterilization methods on the structure and properties of poly(D, L lactic-co-glycolic acid) scaffolds. *Tissue Eng* 2006, **12**:2717–2727.
163. Quadbeck P, Hauser R, Kümmel K, Standke G, Stephani G, Nies B, Rößler S, Wegener B: Iron based cellular metals for degradable synthetic bone replacement. PM2010 World Congress. *PM Biomaterials* 2010, 1–8.
164. Yue S, Lee PD, Pooologasundarampillai G: Evaluation of 3D bioactive glass scaffolds dissolution in a perfusion flow system with X-ray microtomography. *Acta Biomater* 2011, **7**:2637–43.
165. Heinemann C, Heinemann S, Bernhardt A, Worch H, Hanke T: Novel textile chitosan scaffolds promote spreading, proliferation, and differentiation of osteoblasts. *Biomacromolecules* 2008, **9**:2913–2920.
166. Heinemann C, Heinemann S, Lode A, Bernhardt A, Worch H, Hanke T: In vitro evaluation of textile chitosan scaffolds for tissue engineering using human bone marrow stromal cells. *Biomacromolecules* 2009, **10**:1305–1310.
167. Heinemann C, Heinemann S, Bernhardt A, Lode A, Worch H, Hanke T: In vitro osteoclastogenesis on textile chitosan scaffold. *Eur Cell Mater* 2010, **19**:96–106.
168. Rath SN, Strobel LA, Arkudas A, Beier JP, Maier AK, Greil P, Horch R, Kneser U: Osteoinduction and survival of osteoblasts and bone marrow stromal cells in 3D biphasic calcium phosphate scaffolds under static and dynamic culture conditions. *J Cellular Molecular Medicine* 2012, doi:10.1111/j.1582-4934.2012.01545.x.
169. Farack J, Wolf-Brandstetter C, Glorius S, Nies B, Standke G, Quadbeck P, Worch H, Scharnweber D: The effect of perfusion culture on proliferation and differentiation of human mesenchymal cells on biocorrodible bone replacement material. *Materials Sci and Eng* 2011, **176**:1776–1722.
170. Volkmer E, Drosse I, Otto S, Stangelmayer A, Stengale M, Kallukalam BC, et al: Hypoxia in static and dynamic 3D culture systems for tissue engineering of bone. *Tissue Eng Part A* 2008, **14**:1331–1340.
171. Ramsauer M, Kunz J, Krause D, Dermietzel R: Regulation of a blood-brain barrier-specific enzyme expressed by cerebral pericytes (pericytic aminopeptidase N/pAPN) under cell culture conditions. *J Cereb Blood Flow Metab* 1998, **18**:1270–1281.
172. Rüder C, Sauter C, Becker T, Kratz K, Hiebl B, Jung F, Lendlein A, Zohlnhofer D: Viability, proliferation and adhesion of smooth muscle cells and human umbilical vein endothelial cells on electrospun polymer scaffolds. *Microcirc* 2012, **50**:101–12.
173. Blöchl A, Sirrenberg C: Neurotrophins stimulate the release of dopamine from rat mesencephalic neurons via Trk and p75Lnr receptors. *J Biol Chem* 1996, **271**:21100–21107.
174. Canossa M, Griesbeck O, Berninger B, Campana G, Kolbeck R, Thoenen H: Neurotrophin release by neurotrophins: implications for activity-dependent neuronal plasticity. *Proc Natl Acad Sci USA* 1997, **94**:13279–13286.
175. Hoener MC: Role played by sodium in activity-dependent secretion of neurotrophins - revisited. *Eur J Neurosci* 2000, **12**:3096–3106.
176. Griesbeck O, Canossa M, Campana G, Gärtner A, Hoener MC, Nawa H, et al: Are there differences between the secretion characteristics of NGF and BDNF? Implications for the modulatory role of neurotrophins in activity-dependent neuronal plasticity. *Microsc Res Tech* 1999, **45**:262–275.
177. Constantinescu R, Constantinescu AT, Reichmann H, Janetzky B: Neuronal differentiation and long-term culture of the human neuroblastoma line SH-SY5Y. *J Neural Transm Suppl* 2007, **72**:17–28.
178. Yamagishi K, Enomoto T, Ohmiya Y: Perfusion-culture-based secreted bioluminescence reporter assay in living cells. *Anal Biochem* 2006, **345**:15–21.
179. Kalarmarz-Kubiak H, Gozdowska M, Nietrzeba M, Kulcykska E: A novel approach to ATV and IT studies in fish brain and pituitary: in vitro perfusion technique. *J Neuroscience Methods* 2011, **199**:56–61.

doi:10.1186/2001-1326-1-22

Cite this article as: Minuth and Denk: Supportive development of functional tissues for biomedical research using the MINUSHEET® perfusion system. *Clinical and Translational Medicine* 2012 **1**:22.

Submit your manuscript to a SpringerOpen® journal and benefit from:

- Convenient online submission
- Rigorous peer review
- Immediate publication on acceptance
- Open access: articles freely available online
- High visibility within the field
- Retaining the copyright to your article

Submit your next manuscript at ► springeropen.com